

# SCALTEC

## SPB

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## Intended Use

The SPB Series from Sartorius offers precision balances/scales with capacities ranging from 1 mg to 6 kg.

These balances meet the highest requirements for the accuracy and reliability of weighing results through the following features:

- Efficient filtering-out of unfavorable ambient conditions such as vibrations, drafts, etc.
- Stable and reproducible results
- Rugged, durable weighing system

SPB balances save work and speed up simple routine applications through:

- Ultrafast response times
- Easy operation

You can also choose from the following extra functions for simple applications:

- Toggling between weight units
- Net-total formulation (tare memory)
- Counting
- Weighing in percent
- Averaging

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## Warning and Safety Information

Read these operating instructions thoroughly before using your balance/scale to prevent damage to the equipment. Keep these instructions in a safe place.

Follow the instructions below to ensure safe and trouble-free operation of your balance:

- ⚠ Make sure that the voltage rating printed on the AC adapter is identical to your local line voltage
- ⚠ Use only a commercially available 9-volt battery or rechargeable battery
- ⚠ Do not use this balance/scale in a hazardous area/location
- The only way to turn the power off completely is to disconnect the AC adapter
- Connect only SCALTEC accessories and options, as these are optimally designed for use with your balance.
- Protect the AC adapter from contact with liquid.

Do not open the balance housing. If the seal is broken, this will result in forfeiture of all claims under the manufacturer's warranty.

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### Warranty

Do not miss out on the benefits of our extended warranty. Please complete the warranty registration card and return it to the address listed below (see supplement).

In the unlikely event that your balance should require servicing within the 3-year warranty period, please contact your local SCALTEC dealer or:

SCALTEC INSTRUMENTS GmbH  
Rudolf-Diesel-Strasse 1  
37308 Heiligenstadt, Germany  
Phone: (0)3606-6701-0  
Fax: (0)3606-6701-60

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## Storage and Shipping Conditions

Do not expose the balance/scale to shocks, vibrations, moisture or extreme temperatures.

### Unpacking the Balance/Scale

- After unpacking the balance, check it immediately for any visible damage as a result of rough handling during shipment.
- If you see any sign of damage, proceed as directed in the chapter entitled "Care and Maintenance," under the section on "Safety Inspection."

Save the box and all parts of the packaging until you have successfully installed your balance. Only the original packaging provides the best protection for shipment. Before packing your balance, unplug all connected cables to prevent damage.

### Equipment Supplied

The equipment supplied includes the components listed below:

- Balance
- Weighing pan/load plate
- Pan support
- Glass draft shield (only with model SPB42)
- Anti-theft locking device
- AC adapter
- Operating instructions

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
### **Installation Instructions**

Your balance is designed to provide reliable weighing results under normal ambient conditions. When choosing a location to set up your balance, observe the following so that you will be able to work with added speed and accuracy:

- Set up the balance on a stable, even surface
- Avoid placing the balance in close proximity to a heater or otherwise exposing the balance to heat or direct sunlight
- Protect the balance from drafts that come from open windows or doors
- Avoid exposing the balance to extreme vibrations during weighing
- Protect the balance from aggressive chemical vapors
- Do not expose the balance to extreme moisture over long periods

Conditioning the Balance:


Moisture in the air can condense on the surfaces of a cold balance whenever it is brought into a substantially warmer place. If you transfer the balance to a warmer area, make sure to condition it for about 2 hours at room temperature, leaving it unplugged from AC power.

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
## Setting up the Balance

Balances with a Glass Draft Shield:

- Place the components listed below on the balance in the order given:
  - Shield ring – place it on the balance so that the edge for fitting the glass draft shield faces upwards and turn it until it is firmly in place.
  - Pan support                      – Glass draft shield
  - Weighing pan                      – Draft shield cover

Balances with a Round Weighing Pan

- Place the components listed below on the balance in the order given:
  - Pan support
  - Weighing pan/load plate

Balances with a Rectangular Weighing Pan

- Place the weighing pan/load plate on the balance



## Connecting the Balance to AC Power/Safety Precautions

Use only original Sartorius AC adapters:

- for Europe: 6971948    – for the RSA: 6971949
- for the US: 6971947    – for Australia: 6971950
- Insert the right-angle plug into the jack
- The AC adapter rated to Class 2 can be plugged into any wall outlet without requiring any additional safety precautions

The ground is connected to the balance housing, which can be additionally grounded for operation.

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
### Using a (Rechargeable) Battery

(not for models SPB42, SPB52)

- Battery is not included with the equipment supplied
- ⚠ Use only a commercially available (rechargeable) 9-volt battery
- ⚠ When using a rechargeable battery, always use an external recharging device to recharge the battery
- Lay the balance on its side
- Open the battery compartment: Lift the compartment cover
- Insert the 9-volt (rechargeable) battery in the compartment
- Make sure to connect the positive and negative poles correctly
- ⚠ All used batteries are classified as waste that requires special handling (not "household" waste). Dispose of rechargeable batteries in accordance with the applicable special waste disposal regulations.
- Close the battery compartment: Press down on the cover until it clicks into place

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### Leveling the Balance

(only for models SBP42, SBP52)

Level the balance any time you set it up in a new location.

Use only the 2 front feet of the balance for leveling.

- Turn the 2 rear feet until they are in position (only on models SBP52)
  - Turn the 2 front feet as shown here in the illustration until the air bubble is centered in the level indicator
- > In most cases this will require several adjustment steps

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## Basic Weighing Function

### Available Features

- Taring the balance  
You can tare the balance within the entire weighing range.

### Preparation

- Turn on the balance: Press [ON/OFF]
- To change configurations: See the chapter entitled "Configuring the Balance"
- To load factory-set configurations:  
See "Configuring the Balance," parameter 9 – 1
- To tare the balance: Press [TARE]

Additional Functions

- To turn off the balance: Press [ON/OFF]

### Example

Basic weighing

Menu code settings: Factory-set codes

Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
Self-test is performed		
2. Place container on balance scale (here: 52 g)		+ 52.0 g
3. Tare the balance	[TARE]	+ 0.0 g
4. Place sample in container on balance (here: 150.2 g)		+ 150.2 g

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## Calibration/Adjustment\*

Always calibrate/adjust the balance after setting it up in a new location.

### Available Features

Calibration/adjustment can only be performed when

- there is no load on the balance,
- the balance is tared,
- the internal signal is stable.

If these conditions are not met, an error message is displayed. The weight required for calibration/adjustment is displayed (see "Accessories" for calibration weights).

Calibration/adjustment can be performed with any of three weight units:

g, kg, or lb (menu code 1. 4. x)

Factory setting: Grams (menu code 1. 4. 1)

### Example

Calibrate the balance


Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
2. Tare the balance	[TARE]	+ 0.0 g
3. Begin calibration Calibration weight is displayed without weight unit (here: 1000 g)	[CAL/CF]	+ 1000.0 g
4. Place the indicated calibration weight on the balance		+ 150.2 g
After calibration, the calibration weight is displayed with wt. unit		+ 1000.0 g
5. Remove the calibration weight		0.0 g

\* "Calibration" technically means to determine the difference between the balance readout and the actual weight on the pan to determine the accuracy. Adjustment means to bring a balance into the state of accuracy required for its use. Therefore, "calibration," as used in this manual, actually means "adjustment."

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
### Counting (Code 2. 1. 4)

#### Purpose

The counting program allows you to count the number of parts that have approximately the same weight.

#### Available Features

- Store the current weight value to have it loaded as the preset reference sample quantity next time you initialize the Counting application
- The reference sample quantity can be changed in the operating menu: See "Configuring the Balance"
- The average piece weight is automatically output via the optional data interface port after initialization, if the menu code for "Printout with data ID codes" is set
- Press [FUNCTION] to toggle the display between piece count and weight

#### Factory Settings

Reference sample quantity: 10 (menu code 3. 3. 2)

#### Preparation

- Configure the Counting application in the operating menu: See "Configuring the Balance"  
Menu code 2. 1. 4 Counting
  - Reference sample quantity:
    - Menu code 3. 3. 1 5 pcs
    - Menu code 3. 3. 2 10 pcs
    - Menu code 3. 3. 3 20 pcs
    - Menu code 3. 3. 4 50 pcs
    - Menu code 3. 3. 5 100 pcs
- See also "Configuring the Balance"

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## Example

Determine an unknown piece count; weigh the preset reference sample quantity

Settings (changes in the factory settings required for this example):

Menu: Application program: Counting (2. 1. 4)  
Reference sample quantity: 20 pcs (3. 3. 3)

Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
2. Tare the balance	[TARE]	0.0 g
3. Display the reference sample quantity (here: 20 pcs)	[FUNCTION] >2 sec.	rEF 20 (briefly)
4. Place the reference sample quantity (20 pcs) on the balance (here: 66 g)		+ 66.0 g
5. Start the application; if the print format is set to include data ID codes, the following is printed	[FUNCTION]	+ 20 pcs wRef + 3.300 g
6. Weigh uncounted parts (here: 174 pcs)		+ 174 pcs
7. Display weight	[FUNCTION]	+ 574.2 g
8. Display quantity	[FUNCTION]	+ 174 pcs
9. Unload the balance		0 pcs
10. Delete the reference value	[CAL/CF] >2 sec.	
11. Repeat the procedure starting from Step 6, if desired.		

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## Weighing in Percent (Code 2. 1. 5)

### Purpose

This application program allows you to obtain weight readouts in percent which are in proportion to a reference weight.

### Available Features

- Store the current weight value to have it loaded as the preset reference percentage next time you initialize the Weighing in Percent application
- The reference percentage can be changed in the operating menu: See “Configuring the Balance”
- The reference percentage is automatically output via the optional data interface port after initialization, if the menu code for “Printout with data ID codes” is set
- Press [FUNCTION] to toggle the display between percentage and weight

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## Example

Determine an unknown percentage; store the weight on the balance as a reference percentage

Settings (changes in the factory settings required for this example):

Menu: Application program: Weighing in percent (2. 1. 5)

Menu: Reference percentage 100 % (3. 3. 5)

Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
2. Tare the balance	[TARE]	0.0 g
3. Display the reference percentage	[FUNCTION] >2 sec.	rEF 100
4. Place the reference weight for 100% on the balance (here: 222.5 g)		+ 222.5 g
5. Start application; if the print format is set to include data ID codes, the following is printed	[FUNCTION]	+ 100.00 pcs wx% + 222.500 g
6. Place an unknown weight on the balance (here: 322.5 g)		+ 144.94 %
7. Display weight	[FUNCTION]	+ 322.5 g
8. Display percentage	[FUNCTION]	+ 144.94 pcs
9. Unload the balance		0.00 %
10. Delete the reference percentage	[CAL/CF] >2 sec.	
11. Repeat the procedure starting from Step 6, if desired.		

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## Averaging (Code 2. 1. i2)

### Purpose

Use this program to determine weights under unstable ambient conditions. In this program, the balance calculates the weight as the average value from a defined number of individual weighing operations. These weighing operations are also known as "subweighing operations" or "subweighs."

### Available Features

- The measured result displayed is the arithmetic mean shown in the selected weight unit; a triangle under the plus or minus sign indicates that this is a calculated value
- You can set the number of subweighing operations performed in the operating menu: See "Configuring the Balance"
- Press [FUNCTION] for at least 2 sec. to display the pre-set number of subweighing operations
- Press [FUNCTION] to toggle the display between the calculated result and the weight

### Factory Settings

Number of subweighs for averaging: 10 (3. 3. 2)

### Preparation

- Configure the Averaging application in the operating menu: See "Configuring the Balance"  
Menu code 2. 1. i2 Averaging
- Number of subweighs for averaging:

3. 3. 1	5 subweighs
3. 3. 2	10 subweighs
3. 3. 3	20 subweighs
3. 3. 4	50 subweighs
3. 3. 5	100 subweighs

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20

## Example

Determine the weight of a sample in extremely unstable ambient conditions by calculating the average of 10 subweighing operations.

Settings (changes in the factory settings required for this example):

Menu: Application program: Averaging (menu code 2. 1. i2)

Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
2. Tare the balance	[TARE]	0.0 g
3. Display the number of subweighs (here: 10)	[FUNCTION] >2 sec.	10 (briefly)
4. Place sample on the balance (weight readout fluctuates; here: around 275 g)		+ 8888
5. Start measurement	[FUNCTION]	+ 8888 10 9 8 ... 1
If the print format is set to include data ID codes, the following is printed		+ 275.5 g
6. Unload the balance		RES + 275.5 g
7. Delete the result	[CAL/CF] >2 sec.	+ 275.5 g
8. Repeat the procedure starting from Step 4, if desired.		

## Netto-Total/ Zweiter Taraspeicher (Code 2. 1.3)

## Net-Total Formulation/ Second Tare Memory (Code 2. 1.3)

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With this application program you can weigh in components for formulation of a mixture.

### Preparation

Configure the Net-total formulation/Second tare memory application in the operating menu: See "Configuring the Balance" Menu code 2. 1.3

### Example

Step	Key (or instruction)	Display
1. Turn on the balance	[ON/OFF]	
2. Place an empty container on the balance		+ 65.0 g
3. Tare the balance	[TARE]	
4. Add the first component		+ 120.5 g
5. Store the first component weight If the print format is set to include data ID codes, the following is printed	[FUNCTION]	0.0 g <sub>NET</sub> N1 + 120.5 g
6. Add the next component		+ 70.5 g
7. Store the 2nd component weight	[FUNCTION]	0.0 g <sub>NET</sub>
8. Add further components, if desired	As described for Steps 5 and 6	
9. Display total weight and fill to desired final weight	[CAL/CF] >2 sec.	+ 191.0 g

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## Toggle between Weight Units

With this application program you can toggle the display of a weight value back and forth between two weight units.

Configure the "Toggle weight units" application in the operating menu:  
See "Configuring the Balance" menu code 2. 1. 2 Toggle weight units

Menu code		Unit	Conversion factor	Abbr. on printout
1. 7. 2*	3. 1. 2*	Grams	1.00000000000	g
1. 7. 3	3. 1. 3	Kilograms	0.00100000000	kg
1. 7. 4	3. 1. 4	Carast	5.00000000000	ct
1. 7. 5	3. 1. 5	Pounds	0.00220462260	lb
1. 7. 6	3. 1. 6	Ounces	0.03527396200	oz
1. 7. 7	3. 1. 7	Troy ounces	0.03215074700	ozt
1. 7. 8	3. 1. 8	Hong Kong taels	0.02671725000	tlh
1. 7. 9	3. 1. 9	Singapur taels	0.02645544638	tls
1. 7. 10	3. 1. 10	Taiwanese taels	0.02666666000	tlt
1. 7. 11	3. 1. 11	Grains	15.43235835000	GN
1. 7. 12	3. 1. 12	Pennyweights	0.64301493100	dwt
1. 7. 13	3. 1. 13	Milligrams	1000.00000000000	mg
1. 7. 14	3. 1. 14	Parts per pound	1.12876677120	/lb
1. 7. 15	3. 1. 15	Chinese taels	0.02645547175	tlc
1. 7. 16	3. 1. 16	Mommies	0.26670000000	mom
1. 7. 17	3. 1. 17	Austrian carats	5.00000000000	K
1. 7. 18	3. 1. 18	Tola	0.08573333810	tol
1. 7. 19	3. 1. 19	Baht	0.06578947437	bat
1. 7. 20	3. 1. 20	Mesghal	0.21700000000	MS

\* = Factory setting

### Funktion

- To toggle the display between the 1st and 2nd weight units:  
Press the [FUNCTION] key

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Error codes are shown in the main display for approx. 2 seconds, after which the program automatically returns to the weighing mode.

Display	Cause	Solution
No segments appear on the display	No AC power is available The AC adapter is not plugged in Battery is dead	Check the AC power supply Plug in the AC adapter Replace the battery Recharge the battery using an external recharging device
H	The load exceeds the balance capacity	Unload the balance
L	The load plate is not in place Something is touching the load plate	Place the load plate on the balance Move the object that is touching the load plate
E 0 1	Display capacity exceeded: Value to be output cannot be shown on the display	Decrease the weight on the scale
E 0 2	Calibration parameter not met; e.g.: – balance not zeroed – balance is loaded	Calibrate only when zero is displayed Press <b>[TARE]</b> to tare the balance Unload the balance
E 1 0	The <b>[TARE]</b> key is blocked when there is data in the second tare memory (net-total) – only 1 tare function can be used at a time	Press <b>[CAL/CF]</b> >2 sec. to clear the tare memory and release the tare key
E 1 1	Value input is not allowed for second tare memory	Press <b>[TARE]</b>
E 2 2	Weight is too light or there is no sample on the balance	Increase the reference quantity
E 3 0	Interface port for printer output is blocked	Contact your local SCALTEC Service Center
Max. weighing range is less than indicated under "Specifications"	The balance was turned on without the weighing pan in place	Place the weighing pan on the balance and press <b>[ON]</b> to turn on the balance
The weight readout is obviously wrong	The balance has not been calibrated/adjusted The balance was not tared before weighing	Calibrate/adjust the balance Tare before weighing

If any other errors occur, contact your local SCALTEC Service Center!

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### Service

Regular servicing by a SCALTEC technician will extend the service life of your balance and ensure its continued weighing accuracy.

### Repairs

Repair work must be performed by trained service technicians. Any attempt by untrained persons to perform repairs may lead to hazards for the user.

### Cleaning

- ⚠ Disconnect the balance from the AC adapter and unplug any data cables that are connected to the balance
- ⚠ Make sure that no liquid enters the balance housing
- ⚠ Do not use any aggressive cleaning agents (solvents or similar agents)
- Clean the balance using a piece of cloth which has been wet with a mild detergent (soap)
- After cleaning, wipe down the balance with a soft, dry cloth

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### Safety Inspection

If there is any indication that safe operation of the balance with the AC adapter is no longer warranted:

- Turn off the power and disconnect the equipment from AC power immediately
- > Lock the equipment in a secure place to ensure that it cannot be used for the time being

Safe operation of the balance with the AC adapter is no longer ensured when:

- there is visible damage to the AC adapter
- the AC adapter no longer functions properly
- The AC adapter has been stored for a relatively long period under unfavorable conditions

In this case, notify your nearest SCALTEC Service Center. Maintenance and repair work may only be performed by service technicians who are authorized by SCALTEC and who

- have access to the required maintenance manuals
- have attended the relevant service training courses

### Instructions for Recycling the Packaging

To ensure safe shipment, your balance has been packaged using environmentally friendly materials. After successful installation of the balance, you should return this packaging for recycling.

For information on recycling options, including recycling of old weighing equipment, contact your municipal waste disposal center or local recycling depot.

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### Specifications

General specifications:

AC power source/power requirements		AC adapter, 230 or 115 V, +15%...-20%
Frequency		48–60 Hz
Power requirements, direct current	V	10 to 20
Power consumption (average)	W	0.75
Operating temperature range		-10...+40 °C (263...313 K, 32 °F... 104 °F)
Adaptation to ambient conditions		By selection of 1 of 4 optimized filter levels
Display update (depends on the filter level selected)		0.1–0.4
Dimensions (W × D × H)	mm	200 × 270 × 70 (without glass draft shield cylinder)

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Specifications for the individual models:

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### Specifications for the individual models:

Model		SPB62	SPB63	SPB64	SPB71
Weighing capacity	g	3,100	1,500	610	6,000
Readability	g	0.1	0.1	0.1	1
Tare range (subtractive)	g	3,100	1,500	610	6
Repeatability	≤±g	0.1	0.1	0.1	1
Linearity	≤±g	0.2	0.2	0.2	2
Sensitivity drift within +10...+30 °C/+50...+86 °F	≤±/K	1·10 <sup>-5</sup>	2·10 <sup>-5</sup>	5·10 <sup>-5</sup>	5·10 <sup>-5</sup>
Response time (average)	s	2	2	2	2
External calibration weight (of at least accuracy class...)	g lb	1,000 (F2) 2,000 (F2)	1,000 (M1) 2,000 (M1)	500 (M1) 1,000 (M1)	5,000 (M1) 10,000 (M1)
Net weight, approx.	kg	1.4	1.4	1.4	1.4
Pan size	mm	150 Ø			
Hours of operation:					
– with 9-volt battery, approx.	h	11	11	11	11
– with fully charged rechargeable 9-volt battery, approx.	h	2.5	2.5	2.5	2.5
Hours of operation with the SBP05 rechargeable battery pack	h	25			

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The C marking affixed to the equipment indicates that the equipment meets the requirements of the following Directive(s):

### **Council Directive 89/336/EEC "Electromagnetic compatibility (EMC)"**

Applicable European Standards:

Limitation of emissions:

EN 50081-1 Residential, commercial and light industry

EN 50081-2 Industrial environment

Defined immunity to interference:

EN 50082-1 Residential, commercial and light industry

EN 50082-2 Industrial environment

Important Note:

The operator shall be responsible for any modifications to SCALTEC equipment and for any connections of cables or equipment not supplied by SCALTEC and must check and, if necessary, correct these modifications and connections. On request, SCALTEC will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

### **Council Directive 73/23/EEC**

#### **"Electrical equipment designed for use within certain voltage limits"**

Applicable European Standards:

EN 60950 Safety of information technology equipment including electrical business equipment

EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use  
Part 1: General requirements

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.

---

**SCALTEC  
Instruments GmbH**

Rudolf-Diesel-Strasse 1  
37308 Heiligenstadt, Germany  
Phone (3606) 6701-0  
Telefax (3606) 6701-60

Copyright by Scaltec Instruments GmbH, Heiligenstadt, Germany.  
All rights reserved.

No part of this publication may be reprinted or translated in  
any form or by any means without the prior written permission  
of Scaltec Instruments GmbH.

The status of the information, specifications and illustrations  
in this manual is indicated by the date given below.

Scaltec Instruments GmbH reserves the right to make changes to the technology,  
features, specifications, and design of the equipment without notice.

Status: December 1998, Scaltec Instruments GmbH, Heiligenstadt, Germany.

**98648-008-05**

Printed in Germany  
on paper that has been bleached  
without any use of chlorine · W198-A00 · M  
Publication No. WXX6038-a98121

***SCALTEC***